

MyTana camera/locator model repairability status

These are the models of diagnostic equipment we are currently able to repair.



MS11-NG2




MSA-NG2



MSA-NG PC

(no monitor)



DrainTracker



PGR200

200' push-rod in reel and frame



PGR400

400' push-rod in reel and frame



NGX10 or MY30

200' push-rod in reel and frame

MS11-NG1 <- no WiFi -> MSA-NG1

MS11-COM <- analog -> MSA-COM

150' or 200' push-rod

Camera Heads



1 1/2" head



1 1/16" mini head

wired to pushrod or detachable with either:



NG (digital) connector



4-pin (analog) connector

Control Box/ Monitor



Push-rod
If sending just push rod (out of reel) for repair, specify color and length in form's Comments box



CBX17

(for push-rod with NG connection end)

ACB18

(for push-rod with pin connection end)



CBX15

(for push-rod with pin connection end)

Locators



MT 512+ Locator



LineFinder Transmitter

We are NOT able to repair these legacy models.



1 5/8" Camera head



2 1/8" Camera head



ACB13 Control Box



DCB12 Control Box



ML2011 Locator



ACCU-STIC Locator



PT20 Locator